### **MEDICAL GENOMICS LABORATORY**

## ARPKD Known PKHD1 Mutation Testing (KT2)

## **Ordering Information**

#### Acceptable specimen types:

- Fresh blood sample (3-6 ml EDTA; no time limitations associated with receipt)
- Oragene 575 saliva kit (provided by the MGL)
- DNA sample (25ul volume 3ug, O.D. value at 260:280 ≥1.8)

#### Turnaround time:

15 working days

### Price, CPT codes, and Z code:

Sequencing 1 exon: \$250 (USD – institutional/self-pay); CPT: 81403

Sequencing 2 exons: \$340 (USD – institutional/self-pay); CPT: 81403 (x2)

MLPA: \$250 (USD – institutional/self-pay); CPT: 81407

Z code: ZB67K

#### Candidates for this test:

ARPKD testing services are offered to established patients and families with previously identified *PKHD1* variants

### Specimen shipping and handling:

- Please find acceptable specimen type above.
- All submitted specimens must be sent at room temperature. DO NOT ship on ice.
- Specimens must be packaged to prevent breakage and absorbent material must be

included in the package to absorb liquids in the event that breakage occurs. Also, the

**MEDICAL GENOMICS LABORATORY** 

package must be shipped in double watertight containers (e.g. a specimen pouch + the

shipping company's diagnostic envelope).

• To request a sample collection kit, please visit the website or email

medgenomics@uabmc.edu to complete the specimen request form.

• Please contact the MGL (via email at medgenomics@uabmc.edu, or via phone at 205-934-

5562) prior to sample shipment and provide us with the date of shipment and tracking

number of the package so that we can better ensure receipt of the samples.

Required forms:

• Test Requisition Form

• Form for Customs (for international shipments)

Note: Detailed and accurate completion of this document is necessary for reporting

purposes. The Medical Genomics Laboratory issues its clinical reports based on the

demographic data provided by the referring institution on the lab requisition form. It is the

responsibility of the referring institution to provide accurate information. If an amended

report is necessary due to inaccurate or illegible documentation, additional reports will be

drafted with charge.

Requests for testing may not be accepted for the following reasons:

• No label (patients full name and date of collection) on the specimens

No referring physician's or genetic counselor's names and addresses

No billing information

DNA samples must be extracted in a CLIA or equivalent certified lab

For more information, test requisition forms, or sample collection and mailing kits, please call:

205-934-5562.

Phone: (205) 934-5562

**MEDICAL GENOMICS LABORATORY** 

Disorder Background

Autosomal Recessive Polycystic Kidney Disease (ARPKD) is characterized by enlarged cystic

kidneys and hepatic fibrosis. The diagnosis is often made pre- or neonatally, but some patients

are diagnosed later in life. The severity varies widely, with a high mortality rate in the first

months of life. ARPKD is one of the more common hereditary childhood nephropathies with an

estimated incidence of 1:20,000-1:40,000. The carrier frequency in the general population is

estimated to be 1 in 70 to 1 in 100. Variants in PKHD1 are scattered throughout the gene. Most

families carry their own "private" variants. For more information on the condition please refer

to the review on the GeneTests website and Online Mendelian Inheritance in Man.

The gene for ARPKD, PKHD1 (Polycystic Kidney and Hepatic Disease 1), resides on chromosome

6p21-p12, spans 470 kb of genomic DNA and is the only gene known to be associated with the

wide clinical spectrum of autosomal recessive polycystic kidney disease. 86 exons have been

identified and multiple alternative transcripts are known. Over 300 variants have been

reported. Missense, nonsense, frameshift, splicing and multi-exon deletions can occur and are

located throughout the length of the gene, with no major variant hotspots known, as shown in

the PKHD1 variant database.

**Test Description** 

We offer targeted detection of a previously characterized pathogenic variant within the family.

From a fresh EDTA blood sample, DNA is extracted directly and the target region is amplified

and directly sequenced. To offer this testing service, the proband's variant must be identified

by our laboratory before testing relatives.

REFERENCES available on website.

Other related testing options:

ARPKD PKHD1 Prenatal Targeted Testing (PT2)

www.uab.edu/medicine/genetics/medical-genomics-laboratory

# **LAS MEDICAL GENOMICS LABORATORY**

 ARPKD Known PKHD1 Mutation Linkage Analysis for Informativity (PKDL) & Prenatal Linkage (PKDPL)